

Compound (Metabolite, Biomarker)	Matrix	Method	LLOQ	Status*	Photolabile
Abacavir	Plasma	LC-MS/MS	30 ng/mL	Active	No
Abiraterone	Plasma	LC-MS/MS	0.2 ng/mL	Active	No
Acetylcysteine	Plasma	LC-MS/MS	10 ng/mL	Active	No
Aminosalicilic acid	Plasma	LC-MS/MS	50 ng/mL	Active	No
Amlodipine/ Losartan/ Losartan carboxylic acid	Plasma	LC-MS/MS	0.1 ng/mL 1 ng/mL 5 ng/mL	Active	No
Amoxicillin/ Clavulanic acid	Plasma	LC-MS/MS	50 ng/mL 25 ng/mL	Active	No
Aripiprazole	Plasma	LC-MS/MS	0.1 ng/mL	Active	No
Cefixime	Plasma	LC-MS/MS	40 ng/mL	Active	No
Celecoxib	Plasma	LC-MS/MS	10 ng/mL	Active	No
Clofazimine	Plasma	LC-MS/MS	4 ng/mL	Active	No
Clopidogrel	Plasma	LC-MS/MS	0.02 ng/mL	Active	No
Cycloserine	Plasma	LC-MS/MS	50 ng/mL	Active	No
Donepezil	Plasma	LC-MS/MS	0.1 ng/mL	Active	No
Dabigatran	Plasma	LC-MS/MS	0.5 ng/mL	Active	No
Dasatinib	Plasma	LC-MS/MS	2 ng/mL	Active	No
Deferasirox	Plasma	LC-MS/MS	5 ng/mL	Active	No
Dicholine succinate	Plasma	LC-MS/MS	50 ng/mL	Active	No
Dienogest	Plasma	LC-MS/MS	5 ng/mL	Active	No
Diosmin/ Hesperidin	Plasma	LC-MS/MS	0.02 ng/mL 0.1 ng/mL	Active	No
Duloxetine	Plasma	LC-MS/MS	0.5 ng/mL	Active	No
Efavirenz	Plasma	LC-MS/MS	1 ng/mL	Active	No
Emtricitabine	Plasma	LC-MS/MS	10 ng/mL	Active	No
Entecavir	Plasma	LC-MS/MS	0.5 ng/mL	Active	No
Epoetin beta	Plasma	LC-MS/MS	20 ng/mL	Active	No
Erlotinib	Plasma	LC-MS/MS	10 ng/mL	Active	No
Escitalopram	Plasma	LC-MS/MS	0.5 ng/mL	Active	No
Esomeprazole	Plasma	LC-MS/MS	0.1 ng/mL	Active	No
Exemestane	Plasma	LC-MS/MS	0.05 ng/mL	Active	No
Ezetimibe	Plasma	LC-MS/MS	0.25 ng/mL	Active	No
Febuxostat	Plasma	LC-MS/MS	10 ng/mL	Active	No
Fingolimod	Plasma	LC-MS/MS	0.3 ng/mL	Active	No
Fludarabine	Plasma	LC-MS/MS	0.2 ng/mL	Active	No
Gefitinib	Plasma	LC-MS/MS	0.5 ng/mL	Active	No
Hydrochlorothiazide/ Losartan/ Losartan carboxylic acid	Plasma	LC-MS/MS	1 ng/mL 1 ng/mL 5 ng/mL	Active	No
Hydrochlorothiazide/ Telmisartan	Plasma	LC-MS/MS	1 ng/mL 1 ng/mL	Active	No
Ibrutinib	Plasma	LC-MS/MS	5 ng/mL	Active	No
Ibuprofen	Plasma	LC-MS/MS	150 ng/mL	Active	No
Imatinib	Plasma	LC-MS/MS	10 ng/mL	Active	No

Indapamide/ Perindopril/ Perindoprilat	Plasma	LC-MS/MS	0.3 ng/mL 0.4 ng/mL 0.4 ng/mL	Active	No
Lacosamide	Plasma	LC-MS/MS	10 ng/mL	Active	No
Lamivudine	Plasma	LC-MS/MS	25 ng/mL	Active	No
Lamotrigine	Plasma	LC-MS/MS	10 ng/mL	Active	No
Lapatinib	Plasma	LC-MS/MS	5 ng/mL	Active	No
Lenalidomide	Plasma	LC-MS/MS	10 ng/mL	Active	No
Levofloxacin	Plasma	LC-MS/MS	0.4 ng/mL	Active	No
Linezolid	Plasma	LC-MS/MS	50 ng/mL	Active	No
Lisinopril	Plasma	LC-MS/MS	1 ng/mL	Active	No
Lopinavir/ ritonavir	Plasma	LC-MS/MS	62.5 ng/mL 12.5 ng/mL	Active	No
Loratadine	Plasma	LC-MS/MS	0.1 ng/mL	Active	No
Metformin	Plasma	LC-MS/MS	5 ng/mL	Active	No
Metronidazole	Plasma	LC-MS/MS	0.1 ng/mL	Active	Yes
Mexidol (Emoxypine)	Plasma	LC-MS/MS	0.5 ng/mL	Active	No
Montelukast	Plasma	LC-MS/MS	10 ng/mL	Active	No
Mycophenolic acid	Plasma	LC-MS/MS	100 ng/mL	Active	No
Nilotinib	Plasma	LC-MS/MS	2.4 ng/mL	Active	No
Nitisinone	Plasma	LC-MS/MS	20 ng/mL	Active	No
Olanzapine	Plasma	LC-MS/MS	0.1 ng/mL	Active	No
Olmesartan medoxomil	Plasma	LC-MS/MS	5 ng/mL	Active	No
Oseltamivir/ Oseltamivir carboxilate	Plasma	LC-MS/MS	0.5 ng/mL 1 ng/mL	Active	No
Oxcarbazepine	Plasma	LC-MS/MS	20 ng/mL	Active	No
Pazopanib	Plasma	LC-MS/MS	5 ng/mL	Active	No
Pregabalin	Plasma	LC-MS/MS	100 ng/mL	Active	No
Quetiapine	Plasma	LC-MS/MS	1 ng/mL	Active	No
Raltegravir	Plasma	LC-MS/MS	5 ng/mL	Active	No
Ritonavir	Plasma	LC-MS/MS	12.5 ng/mL	Active	No
Rivaroxaban	Plasma	LC-MS/MS	1 ng/mL	Active	No
Rosuvastatin	Plasma	LC-MS/MS	0.1 ng/mL	Active	No
Saquinavir	Plasma	LC-MS/MS	5 ng/mL	Active	No
Sertraline	Plasma	LC-MS/MS	0.1 ng/mL	Active	No
Sildenafil	Plasma	LC-MS/MS	1 ng/mL	Active	No
Sorafenib	Plasma	LC-MS/MS	50 ng/mL	Active	No
Sofosbuvir	Plasma	LC-MS/MS	0.25 ng/mL	Active	No
Sparfloxacin	Plasma	LC-MS/MS	5 ng/mL	Active	No
Sunitinib	Plasma	LC-MS/MS	0.2 ng/mL	Active	No
Tamsulosin	Plasma	LC-MS/MS	0.01 ng/mL	Active	No
Tenofovir	Plasma	LC-MS/MS	10 ng/mL	Active	No
Terizidone	Plasma	LC-MS/MS	50 ng/mL	Active	No
Thioctic acid (Lipoic acid)	Plasma	LC-MS/MS	5 ng/mL	Active	No
Umifenovir	Plasma	LC-MS/MS	1 ng/mL	Active	Yes
Valaciclovir	Plasma	LC-MS/MS	2 ng/mL	Active	No

Valganciclovir	Plasma	LC-MS/MS	5 ng/mL	Active	No
Valsartan	Plasma	LC-MS/MS	10 ng/mL	Active	No
Voriconazole	Plasma	LC-MS/MS	100 ng/mL	Active	No
Zidovudine	Plasma	LC-MS/MS	1 ng/mL	Active	No

* Inactive methods and older methods may require a revalidation fee to bring them up-to-date with the current EMA/FDA Guideline